

Coronado sleep chair

Operating instructions



1. Extending the sleep surface:

- Locate the pull strap on the front of the sleep chair.
- Grasp the pull strap firmly.
- Pull the strap steadily and firmly to roll out the first portion of the sleep surface.
- Continue pulling until the first portion is fully extended.



2. Elevating the ottoman:

- Once the first portion is rolled out, locate the second strap on the front interior of the sleep chair.
- Lift up the second strap and pull it forward to elevate the ottoman, making it level with the first portion.



3. Adjusting the headrest:

- To form the headrest, allow the back of the sleep chair to naturally fold over.
- Gently guide and position this portion to be level with the rest of the unit.



4. Making the bed:

- Once the sleep chair is fully extended and the headrest is adjusted, prepare the bed by adding sheets and a pillow.
- Arrange pillows at the headrest to provide additional comfort and support.

NOTE:

Please insure that seat is completely closed before using as a chair. To test place hand on seat and try to slide forward. If seat moves forward push back until fully closed.