

**Patient recliner stationary arms 9300-R**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.25 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Outside back</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	28.25"	<b>Weight</b>	182 lbs.
<b>Base Clearance</b>	2.25"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Patient recliner stationary arms 9300-R-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Arms</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	28.25"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	2.25"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	182 lbs.
		<b>Weight Rating</b>	500 lbs.



**Patient recliner stationary arms and trendelenburg 9317**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.25 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Outside back</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	28.25"	<b>Weight</b>	182 lbs.
<b>Base Clearance</b>	2.25"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Patient recliner stationary arms and trendelenburg 9317-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Arms</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	28.25"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	2.25"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	182 lbs.
		<b>Weight Rating</b>	500 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Standard with stationary, non-pivoting arms</li> <li>Durable steel mechanism (10 year mechanism warranty)</li> <li>Wall-saver design</li> <li>Foot ottoman and back operate independently of one another</li> <li>Dual backrest activation levers</li> <li>Ergonomic arm design provides long-term comfort</li> <li>Dual-sided central locking caster system</li> <li>Durable urethane push handle</li> <li>Foot ottoman handle available on left or right side</li> <li>Perimeter base trim protects upholstery</li> <li>Due to Lasata's design the outside back has limitations in textile applications</li> <li><b>NOTE:</b> Recliner is not intended for patient transport</li> </ul>	<p><b>Functionality</b></p> <ul style="list-style-type: none"> <li>Left or right side facing foot ottoman handle</li> <li>Perimeter base trim (Light Grey - Frost, Dark Grey - Shadow)</li> </ul>
	<p><b>Additional comfort</b></p> <ul style="list-style-type: none"> <li>Heated seat and back (standard or fully retractable cord)</li> <li>Head pillow</li> </ul>
	<p><b>Accessories and protection</b></p> <ul style="list-style-type: none"> <li>Foley hook (Light Grey - Frost, Dark Grey - Shadow)</li> <li>Removable seat cover</li> <li>Moisture barrier</li> </ul>

**Patient recliner pivoting arms 9400-R**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.25 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Outside back</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	28.25"	<b>Weight</b>	182 lbs.
<b>Base Clearance</b>	2.25"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Patient recliner pivoting arms 9400-R-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Arms</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	28.25"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	2.25"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	182 lbs.
		<b>Weight Rating</b>	500 lbs.



**Patient recliner pivoting arms and trendelenburg 9417**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.25 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Outside back</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	28.25"	<b>Weight</b>	182 lbs.
<b>Base Clearance</b>	2.25"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Patient recliner pivoting arms and trendelenburg 9417-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 49.75"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 20"h	<b>COM - Arms</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	28.25"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	2.25"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	182 lbs.
		<b>Weight Rating</b>	500 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Standard with dual-locking pivoting arms</li> <li>Durable steel mechanism (10 year mechanism warranty)</li> <li>Wall-saver design</li> <li>Foot ottoman and back operate independently of one another</li> <li>Dual backrest activation levers</li> <li>Ergonomic arm design provides long-term comfort</li> <li>Dual-sided central locking caster system</li> <li>Durable urethane push handle</li> <li>Foot ottoman handle available on left or right side</li> <li>Perimeter base trim protects upholstery</li> <li>Due to Lasata's design the outside back has limitations in textile applications</li> <li><b>NOTE:</b> Recliner is not intended for patient transport</li> </ul>	<p><b>Functionality</b></p> <ul style="list-style-type: none"> <li>Left or right side facing foot ottoman handle</li> <li>Perimeter base trim (Light Grey - Frost, Dark Grey - Shadow)</li> </ul>
	<p><b>Additional comfort</b></p> <ul style="list-style-type: none"> <li>Heated seat and back (standard or fully retractable cord)</li> <li>Head pillow</li> </ul>
	<p><b>Accessories and protection</b></p> <ul style="list-style-type: none"> <li>Foley hook (Light Grey - Frost, Dark Grey - Shadow)</li> <li>Removable seat cover</li> <li>Moisture barrier</li> </ul>

**Ortho recliner stationary arms 9310-R-OR**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.75 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Outside back</b>	2.25 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	30"	<b>Weight</b>	182 lbs.
<b>Base Clearance</b>	5.125"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Ortho recliner stationary arms 9310-R-OR-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Arms</b>	2.5 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	30"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	5.125"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	182 lbs.
		<b>Weight Rating</b>	500 lbs.



**Ortho recliner stationary arms and trendelenburg 9312-R-OR**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.75 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Outside back</b>	2.25 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	30"	<b>Weight</b>	185 lbs.
<b>Base Clearance</b>	4.75"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Ortho recliner stationary arms and trendelenburg 9312-R-OR-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Arms</b>	2.5 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	30"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	4.75"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	185 lbs.
		<b>Weight Rating</b>	500 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Standard with stationary, non-pivoting arms</li> <li>Durable steel mechanism (10 year mechanism warranty)</li> <li>Wall-saver design</li> <li>Foot ottoman and back operate independently of one another</li> <li>Dual backrest activation levers</li> <li>Ergonomic arm design provides long-term comfort</li> <li>Dual-sided central locking caster system</li> <li>Durable urethane push handle</li> <li>Foot ottoman handle available on left or right side</li> <li>Perimeter base trim protects upholstery</li> <li>Due to Lasata's design the outside back has limitations in textile applications</li> <li><b>NOTE:</b> Recliner is not intended for patient transport</li> </ul>	<p><b>Functionality</b></p> <ul style="list-style-type: none"> <li>Left or right side facing foot ottoman handle</li> <li>Perimeter base trim (Light Grey - Frost, Dark Grey - Shadow)</li> </ul>
	<p><b>Additional comfort</b></p> <ul style="list-style-type: none"> <li>Heated seat and back (standard or fully retractable cord)</li> <li>Head pillow</li> </ul>
	<p><b>Accessories and protection</b></p> <ul style="list-style-type: none"> <li>Foley hook (Light Grey - Frost, Dark Grey - Shadow)</li> <li>Removable seat cover</li> <li>Moisture barrier</li> </ul>

**Ortho recliner pivoting arms 9410-R-OR**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.75 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Outside back</b>	2.25 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	30"	<b>Weight</b>	182 lbs.
<b>Base Clearance</b>	5.125"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Ortho recliner pivoting arms 9410-R-OR-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Arms</b>	2.5 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	30"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	5.125"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	182 lbs.
		<b>Weight Rating</b>	500 lbs.



**Ortho recliner pivoting arms and trendelenburg 9412-R-OR**  
Single textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	4.75 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Outside back</b>	2.25 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	97.5
<b>Top of Arm to Floor</b>	30"	<b>Weight</b>	185 lbs.
<b>Base Clearance</b>	4.75"	<b>Weight Rating</b>	500 lbs.
<b>Extended Length</b>	77.75"		



**Ortho recliner pivoting arms and trendelenburg 9412-R-OR-MT**  
Multiple textiles

<b>Overall Dimensions</b>	32.75"w x 34.5"d x 53"h	<b>COM - Seat and Inside Back</b>	2 yds.
<b>Seat Dimensions</b>	22.75"w x 19.75"d x 22.75"h	<b>COM - Arms</b>	2.5 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1 yds.
<b>Top of Arm to Floor</b>	30"	<b>COM - Outside Back</b>	2 yds.
<b>Base Clearance</b>	4.75"	<b>Cubic Feet</b>	97.5
<b>Extended Length</b>	77.75"	<b>Weight</b>	185 lbs.
		<b>Weight Rating</b>	500 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Standard with dual-locking pivoting arms</li> <li>Durable steel mechanism (10 year mechanism warranty)</li> <li>Wall-saver design</li> <li>Foot ottoman and back operate independently of one another</li> <li>Dual backrest activation levers</li> <li>Ergonomic arm design provides long-term comfort</li> <li>Dual-sided central locking caster system</li> <li>Durable urethane push handle</li> <li>Foot ottoman handle available on left or right side</li> <li>Perimeter base trim protects upholstery</li> <li>Due to Lasata's design the outside back has limitations in textile applications</li> <li><b>NOTE:</b> Recliner is not intended for patient transport</li> </ul>	<p><b>Functionality</b></p> <ul style="list-style-type: none"> <li>Left or right side facing foot ottoman handle</li> <li>Perimeter base trim (Light Grey - Frost, Dark Grey - Shadow)</li> </ul>
	<p><b>Additional comfort</b></p> <ul style="list-style-type: none"> <li>Heated seat and back (standard or fully retractable cord)</li> <li>Head pillow</li> </ul>
	<p><b>Accessories and protection</b></p> <ul style="list-style-type: none"> <li>Foley hook (Light Grey - Frost, Dark Grey - Shadow)</li> <li>Removable seat cover</li> <li>Moisture barrier</li> </ul>

**Bariatric recliner stationary arms 9301-R-OB**  
Single textiles

<b>Overall Dimensions</b>	39.75"w x 34.5"d x 49.75"h	<b>COM - Seat, Inside Back, Arms, Foot Ottoman</b>	5 yds.
<b>Seat Dimensions</b>	29.75"w x 19.75"d x 20"h	<b>COM - Outside back</b>	2.75 yds.
<b>Seat to Top of Back</b>	30"	<b>Cubic Feet</b>	117
<b>Top of Arm to Floor</b>	28.25"	<b>Weight</b>	194 lbs.
<b>Base Clearance</b>	2.25"	<b>Weight Rating</b>	750 lbs.
<b>Extended Length</b>	77.75"		



**Bariatric recliner stationary arms 9301-R-OB-MT**  
Multiple textiles

<b>Overall Dimensions</b>	39.75"w x 34.5"d x 49.75"h	<b>COM - Seat and Inside Back</b>	2.25 yds.
<b>Seat Dimensions</b>	29.75"w x 19.75"d x 20"h	<b>COM - Arms</b>	2 yds.
<b>Seat to Top of Back</b>	30"	<b>COM - Foot ottoman</b>	1.25 yds.
<b>Top of Arm to Floor</b>	28.25"	<b>COM - Outside Back</b>	2.75 yds.
<b>Base Clearance</b>	2.25"	<b>Cubic Feet</b>	117
<b>Extended Length</b>	77.75"	<b>Weight</b>	194 lbs.
		<b>Weight Rating</b>	750 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Rated to 750 lbs.</li> <li>Standard with stationary, non-pivoting arms</li> <li>Durable steel mechanism (10 year mechanism warranty)</li> <li>Wall-saver design</li> <li>Foot ottoman and back operate independently of one another</li> <li>Dual backrest activation levers</li> <li>Ergonomic arm design provides long-term comfort</li> <li>Dual-sided central locking caster system</li> <li>Durable urethane push handle</li> <li>Foot ottoman handle available on left or right side</li> <li>Perimeter base trim protects upholstery</li> <li>Due to Lasata's design the outside back has limitations in textile applications</li> <li><b>NOTE:</b> Recliner is not intended for patient transport</li> </ul>	<p><b>Functionality</b></p> <ul style="list-style-type: none"> <li>Left or right side facing foot ottoman handle</li> <li>Perimeter base trim (Light Grey - Frost, Dark Grey - Shadow)</li> </ul>
	<p><b>Additional comfort</b></p> <ul style="list-style-type: none"> <li>Heated seat and back (standard or fully retractable cord)</li> <li>Head pillow</li> </ul>
	<p><b>Accessories and protection</b></p> <ul style="list-style-type: none"> <li>Foley hook (Light Grey - Frost, Dark Grey - Shadow)</li> <li>Removable seat cover</li> <li>Moisture barrier</li> </ul>

**Standard Textile & Finish Selection**

---

<http://carolina.ofs.com/resources/textiles>  
<http://carolina.ofs.com/resources/finishes>

**Frame Construction**

---

- |            |   |
|------------|---|
| Seat       | • Plywood                                       |
| Back       | • Contoured hardwood frame construction         |
| Lower Arms | • Plywood/hardwood frame construction           |
| Assembly   | • Bolted connections and general grade hardware |

**Foam**

---

- |           |  |
|-----------|--|
| Seat      | • All foam is premium grade CAL 117 certified fire retardant |
| Back      | • All foam is premium grade CAL 117 certified fire retardant |
| Lower Arm | • All foam is premium grade CAL 117 certified fire retardant |

**Suspension**

---

- |      |                                |
|------|--------------------------------|
| Seat | • Webrider pocket coil springs |
| Back | • Synet elastic sheet webbing  |

**Dual-sided Central Locking Casters**

---

- When the lock is engaged, by design, only the right front-facing, and left rear-facing casters lock, securing the recliner in a locked position.
- The left front-facing, and right rear-facing casters remained unlocked.

**Notes and disclaimers**

---

- Recliner is not to be used for patient transport