

Patient chair stationary back 1478
Single textiles

Overall Dimensions	25" w x 26.5" d x 46.5" h	COM 54"	2 yds.
Seat Dimensions	21.5" w x 20" d x 18" h	Cubic Feet	24
Seat to Top of Back	28.5"	Weight	55
Top of Arm to Floor	25.5"	Weight Rating	300 lbs.
Exposed Metal	Luster Grey		



Patient chair stationary back 1478-MT
Multiple textiles

Overall Dimensions	25" w x 26.5" d x 46.5" h	COM 54" - Seat	1 yd
Seat Dimensions	21.5" w x 20" d x 18" h	COM 54" - Back	1.75 yds.
Seat to Top of Back	28.5"	Cubic Feet	24
Top of Arm to Floor	25.5"	Weight	55
Exposed Metal	Luster Grey	Weight Rating	300 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> Stationary back Upholstered seat and back Wall-saver design Clean-thru design Tight cushion Luster Grey powder coat on all exposed metal Black poly arm caps Non-marring adjustable glides Field replaceable components 	Arm caps <ul style="list-style-type: none"> Poly (Black, Taupe, and Grey) Wood (Standard finish offering) Solid surface (Cameo White, Modern White, Medea, and Polaris)
	Frame <ul style="list-style-type: none"> Powder-coated (standard or super matte powder coats)
	Shift rear caster <ul style="list-style-type: none"> Caster body will match the specified frame finish <p>Note: Rear casters are for repositioning only, and are not designed for patient transport. Do not move the chair while occupied.</p>

Patient chair motion back 1488
Single textiles

Overall Dimensions	25" w x 27.5" d x 46.5" h	COM 54"	2 yds.
Seat Dimensions	21.5" w x 20" d x 18" h	Cubic Feet	24
Seat to Top of Back	28.5"	Weight	65
Top of Arm to Floor	25.5"	Weight Rating	300 lbs.
Exposed Metal	Luster Grey		



Patient chair motion back 1488-MOTION-MT
Multiple textiles

Overall Dimensions	25" w x 27.5" d x 46.5" h	COM 54" - Seat	1 yd
Seat Dimensions	21.5" w x 20" d x 18" h	COM 54" - Back	1.75 yds.
Seat to Top of Back	28.5"	Cubic Feet	24
Top of Arm to Floor	25.5"	Weight	65
Exposed Metal	Luster Grey	Weight Rating	300 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> Motion back Upholstered seat and back Wall-saver design Clean-thru design Tight cushion Luster Grey powder coat on all exposed metal Black poly arm caps Non-marring adjustable glides Field replaceable components 	Arm caps <ul style="list-style-type: none"> Poly (Black, Taupe, and Grey) Wood (Standard finish offering) Solid surface (Cameo White, Modern White, Medea, and Polaris)
	Frame <ul style="list-style-type: none"> Powder-coated (standard or super matte powder coats)
	Shift rear caster <ul style="list-style-type: none"> Caster body will match the specified frame finish <p>Note: Rear casters are for repositioning only, and are not designed for patient transport. Do not move the chair while occupied.</p>

Easy access chair 1474
Single textiles

Overall Dimensions	25" w x 27.5" d x 40.5" h	COM 54"	1.25 yds.
Seat Dimensions	21.5" w x 20" d x 23.75" h	Cubic Feet	24
Seat to Top of Back	15.25"	Weight	60
Top of Arm to Floor	31.25"	Weight Rating	500 lbs.
Exposed Metal	Luster Grey		



Easy access chair 1474-MOTION-MT
Multiple textiles

Overall Dimensions	25" w x 27.5" d x 40.5" h	COM 54" - Seat	1 yd
Seat Dimensions	21.5" w x 20" d x 23.75" h	COM 54" - Back	1.75 yds.
Seat to Top of Back	15.25"	Cubic Feet	24
Top of Arm to Floor	31.25"	Weight	60
Exposed Metal	Luster Grey	Weight Rating	500 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> Upholstered seat and back Wall-saver design Clean-thru design Tight cushion Luster Grey powder coat on all exposed metal Black poly arm caps Polyurethane coated steel foot rest Non-marring adjustable glides Field replaceable components 	Arm caps <ul style="list-style-type: none"> Poly (Black, Taupe, and Grey) Wood (Standard finish offering) Solid surface (Cameo White, Modern White, Medea, and Polaris)
	Frame <ul style="list-style-type: none"> Powder-coated (standard or super matte powder coats)
	Foot rest <ul style="list-style-type: none"> Poly (Black, Taupe, and Grey)
	Shift rear caster <ul style="list-style-type: none"> Caster body will match the specified frame finish <p>Note: Rear casters are for repositioning only, and are not designed for patient transport. Do not move the chair while occupied.</p>

Standard Textile & Finish Selection

<http://carolina.ofs.com/resources/textiles>
<http://carolina.ofs.com/resources/finishes>

Frame Construction

Seat	<ul style="list-style-type: none"> Plywood seat EZ frame construction
Back	<ul style="list-style-type: none"> Contoured plywood shape 11-ply
Arms	<ul style="list-style-type: none"> 1" 16 gauge C1018 steel posts 3/4" x 1/2" 12 gauge C1018 rails
Assembly	<ul style="list-style-type: none"> All welded joints are full penetration mig/tig welds All bolted connections are 1/4" 20 grade 5 hardware

Foam

Seat/Back	<ul style="list-style-type: none"> All foam is premium grade CAL 117 certified fire retardant
-----------	--

Suspension

- LEXUS webbing
- Polypropylene/latex/rubber
- 2430 propex (applied to backs)
- 100% polypropylene
- Latex coating is paranol acrylic 763
- Applied with 5/16 stainless steel staples