

**Patient chair stationary back 1478**  
Single textiles

<b>Overall Dimensions</b>	25" w x 26.5" d x 46.5" h	<b>COM 54"</b>	2 yds.
<b>Seat Dimensions</b>	21.5" w x 20" d x 18" h	<b>Cubic Feet</b>	24
<b>Seat to Top of Back</b>	28.5"	<b>Weight</b>	55
<b>Top of Arm to Floor</b>	25.5"	<b>Weight Rating</b>	300 lbs.
<b>Exposed Metal</b>	Luster Grey		



**Patient chair stationary back 1478-MT**  
Multiple textiles

<b>Overall Dimensions</b>	25" w x 26.5" d x 46.5" h	<b>COM 54" - Seat</b>	1 yd
<b>Seat Dimensions</b>	21.5" w x 20" d x 18" h	<b>COM 54" - Back</b>	1.75 yds.
<b>Seat to Top of Back</b>	28.5"	<b>Cubic Feet</b>	24
<b>Top of Arm to Floor</b>	25.5"	<b>Weight</b>	55
<b>Exposed Metal</b>	Luster Grey	<b>Weight Rating</b>	300 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Stationary back</li> <li>Upholstered seat and back</li> <li>Wall-saver design</li> <li>Clean-thru design</li> <li>Tight cushion</li> <li>Luster Grey powder coat on all exposed metal</li> <li>Black poly arm caps</li> <li>Non-marring adjustable glides</li> <li>Field replaceable components</li> </ul>	<b>Arm caps</b> <ul style="list-style-type: none"> <li>Poly (Black, Taupe, and Grey)</li> <li>Wood (Standard finish offering)</li> <li>Solid surface (Cameo White, Modern White, Medea, and Polaris)</li> </ul>
	<b>Frame</b> <ul style="list-style-type: none"> <li>Powder-coated (standard or super matte powder coats)</li> </ul>
	<b>Shift rear caster</b> <ul style="list-style-type: none"> <li>Caster body will match the specified frame finish</li> </ul> <p><b>Note:</b> Rear casters are for repositioning only, and are not designed for patient transport. Do not move the chair while occupied.</p>

**Patient chair motion back 1488**  
Single textiles

<b>Overall Dimensions</b>	25" w x 27.5" d x 46.5" h	<b>COM 54"</b>	2 yds.
<b>Seat Dimensions</b>	21.5" w x 20" d x 18" h	<b>Cubic Feet</b>	24
<b>Seat to Top of Back</b>	28.5"	<b>Weight</b>	65
<b>Top of Arm to Floor</b>	25.5"	<b>Weight Rating</b>	300 lbs.
<b>Exposed Metal</b>	Luster Grey		



**Patient chair motion back 1488-MOTION-MT**  
Multiple textiles

<b>Overall Dimensions</b>	25" w x 27.5" d x 46.5" h	<b>COM 54" - Seat</b>	1 yd
<b>Seat Dimensions</b>	21.5" w x 20" d x 18" h	<b>COM 54" - Back</b>	1.75 yds.
<b>Seat to Top of Back</b>	28.5"	<b>Cubic Feet</b>	24
<b>Top of Arm to Floor</b>	25.5"	<b>Weight</b>	65
<b>Exposed Metal</b>	Luster Grey	<b>Weight Rating</b>	300 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Motion back</li> <li>Upholstered seat and back</li> <li>Wall-saver design</li> <li>Clean-thru design</li> <li>Tight cushion</li> <li>Luster Grey powder coat on all exposed metal</li> <li>Black poly arm caps</li> <li>Non-marring adjustable glides</li> <li>Field replaceable components</li> </ul>	<b>Arm caps</b> <ul style="list-style-type: none"> <li>Poly (Black, Taupe, and Grey)</li> <li>Wood (Standard finish offering)</li> <li>Solid surface (Cameo White, Modern White, Medea, and Polaris)</li> </ul>
	<b>Frame</b> <ul style="list-style-type: none"> <li>Powder-coated (standard or super matte powder coats)</li> </ul>
	<b>Shift rear caster</b> <ul style="list-style-type: none"> <li>Caster body will match the specified frame finish</li> </ul> <p><b>Note:</b> Rear casters are for repositioning only, and are not designed for patient transport. Do not move the chair while occupied.</p>

**Easy access chair 1474**  
Single textiles

<b>Overall Dimensions</b>	25" w x 27.5" d x 40.5" h	<b>COM 54"</b>	1.25 yds.
<b>Seat Dimensions</b>	21.5" w x 20" d x 23.75" h	<b>Cubic Feet</b>	24
<b>Seat to Top of Back</b>	15.25"	<b>Weight</b>	60
<b>Top of Arm to Floor</b>	31.25"	<b>Weight Rating</b>	300 lbs.
<b>Exposed Metal</b>	Luster Grey		



**Easy access chair 1474-MOTION-MT**  
Multiple textiles

<b>Overall Dimensions</b>	25" w x 27.5" d x 40.5" h	<b>COM 54" - Seat</b>	1 yd
<b>Seat Dimensions</b>	21.5" w x 20" d x 23.75" h	<b>COM 54" - Back</b>	1.75 yds.
<b>Seat to Top of Back</b>	15.25"	<b>Cubic Feet</b>	24
<b>Top of Arm to Floor</b>	31.25"	<b>Weight</b>	60
<b>Exposed Metal</b>	Luster Grey	<b>Weight Rating</b>	300 lbs.



Standard Features	Options Available
<ul style="list-style-type: none"> <li>Upholstered seat and back</li> <li>Wall-saver design</li> <li>Clean-thru design</li> <li>Tight cushion</li> <li>Luster Grey powder coat on all exposed metal</li> <li>Black poly arm caps</li> <li>Polyurethane coated steel foot rest</li> <li>Non-marring adjustable glides</li> <li>Field replaceable components</li> </ul>	<b>Arm caps</b> <ul style="list-style-type: none"> <li>Poly (Black, Taupe, and Grey)</li> <li>Wood (Standard finish offering)</li> <li>Solid surface (Cameo White, Modern White, Medea, and Polaris)</li> </ul>
	<b>Frame</b> <ul style="list-style-type: none"> <li>Powder-coated (standard or super matte powder coats)</li> </ul>
	<b>Foot rest</b> <ul style="list-style-type: none"> <li>Poly (Black, Taupe, and Grey)</li> </ul>
	<b>Shift rear caster</b> <ul style="list-style-type: none"> <li>Caster body will match the specified frame finish</li> </ul> <p><b>Note:</b> Rear casters are for repositioning only, and are not designed for patient transport. Do not move the chair while occupied.</p>

**Standard Textile & Finish Selection**

<http://carolina.ofs.com/resources/textiles>  
<http://carolina.ofs.com/resources/finishes>

**Frame Construction**

Seat	<ul style="list-style-type: none"> <li>Plywood seat EZ frame construction</li> </ul>
Back	<ul style="list-style-type: none"> <li>Contoured plywood shape 11-ply</li> </ul>
Arms	<ul style="list-style-type: none"> <li>1" 16 gauge C1018 steel posts</li> <li>3/4" x 1/2" 12 gauge C1018 rails</li> </ul>
Assembly	<ul style="list-style-type: none"> <li>All welded joints are full penetration mig/tig welds</li> <li>All bolted connections are 1/4" 20 grade 5 hardware</li> </ul>

**Foam**

Seat/Back	<ul style="list-style-type: none"> <li>All foam is premium grade CAL 117 certified fire retardant</li> </ul>
-----------	--

**Suspension**

- LEXUS webbing
- Polypropylene/latex/rubber
- 2430 propex (applied to backs)
- 100% polypropylene
- Latex coating is paranol acrylic 763
- Applied with 5/16 stainless steel staples